


ALTERNATING PRESSURE MATTRESS WITH PUMP



T08+P06C

ISO13485 FDA CE MDR MDMA



Schematic diagram of A/B cycle alternating mode

- ▼ Stage IV Intensive care
- ▼ 2-way stretch PU cover with snap fixture, waterproof, fire retardant(option)
- ▼ CPR for quick air release in emergency
- ▼ Ventilated bladder
- ▼ Warranty: 1 year




MATTRESS T08

Size:	80"L*32"~36"W*9"H(20cells)
Bariatric size:	42"/48"(width)
Mattress material:	Nylon TPU or Pure TPU or Nylon PVC
Cover material:	PU
Bottom material:	Dacron PVC
Net weight:	8.35Kg
Weight capacity:	210Kg
*Three layer mattress;	
Alternating mode: A-B or A-B-C or A-B-C~S or A-B~S;	
Pillow function	

PUMP P06C

Power:	110V/60Hz or 220V/50Hz
Air output:	7/12 L/min
Cycle time:	6/9/12/25 mins
Size:	345*115*185mm(L*W*H)
Net Weight:	2.55Kg
*Touch screen	



Option:
Top layer and second layer A-B or A-B-C line alternate, bottom layer and pillow keep static.



The vents on a cyclic alternating air mattress are primarily used to regulate air pressure and ventilation, preventing overheating or moisture buildup.

P06B



Pressure display
Size: 345*115*185mm(L*W*H)
Net Weight: 2.55Kg